

Systematic Fractionation Using Complementary Developing Solvent Technique in Flash Chromatography: A Key Step in Dereplication Strategies for Natural Products

Jason Fauquet¹, Adam Ellatiff¹, Claudio Palmieri¹, Pierre Duez¹, Bertrand Blankert², Mathilde Wells² and Amandine Nachtergaele¹

¹Laboratory of Therapeutic Chemistry and Pharmacognosy, Research Institute for Health Sciences and Technology, Research Institute for Bioscience, Faculty of Medicine, Pharmacy and Biomedical Sciences, University of Mons (UMONS), Belgium

²Laboratory of Pharmaceutical Analysis, Research Institute for Health Sciences and Technology, Research Institute for Bioscience, Faculty of Medicine, Pharmacy and Biomedical Sciences, University of Mons (UMONS), Belgium

What?

COMPLEMENTARY DEVELOPING SOLVENT (CDS)

A trio of 3 solvent mixtures: Low, Medium, and High-Polarity Developing Solvent (LPDS, MPDS, HPDS).

Mixture name	Solvents	Proportion (v/v)
HPDS	Ethanol, DCM, Water, FA	43:43:11:3
MPDS	MTBE, THF, Water, FA	61:36:1.5:1.5
LPDS	Toluene, Ethyl acetate	9:1

- Significantly enhance chromatographic resolution.
- Covering a broad range of selectivities.
- Three complementary chromatographic fingerprints


Originally, developed for **untargeted analysis** in High-Performance Thin-Layer Chromatography (HPTLC).

OUR OBJECTIVE ?

Adapting CDS from HPTLC to Flash chromatography



How?




Using the same System suitability test (SST) covering a wide polarity range for direct comparisons :

Universal HPTLC Mix (UHM)

Label	Name	Log P	HPTLC concentration (mg/mL)	Flash concentration (mg/mL)
A	Guanosine	-1.90	0.5	1
B	Sulisobenzene	0.88*	1	1
C	Thymidine	-0.93	1	1
D	Paracetamol	0.46	1	1
D	Phthalimide	1.15	2	1
F	9-hydroxyfluorene	2.84*	1	1
G	Thioxanthene-9-one	3.90*	0.01	1
H	Octrizole	5.88*	1	1

* Computed values (DrugBank, PubChem or SciFinder)


THE FLASH METHOD WE PROPOSE



Key features


- Same elution order as for HPTLC & improved resolution.
- UHM ensures reproducibility in solvent gradient control & maintains consistent elution patterns
- The applicability of HPDS and LPDS gradients to diverse plant matrices allows systematization of untargeted studies

Note : Unfortunately, MPDS was not successfully developed because of non-elution of compounds



THE WAY TO INCLUDE IT IN YOUR WORKFLOW


(An example with HPDS)



BONUS

A new systematic fractionation strategy !

- Seven isocratic steps by gradients
- One fraction collected by isocratic step - F1 to F7
- Each fraction is enriched with compounds of similar polarity
- Reduced fraction numbers : 1 extract = max 14 fractions VS many more with peak-to-peak collection strategy
- Accelerates bioguided fractionation and untargeted metabolomics workflows



Mass Spectrometry (MS) Nuclear Magnetic Resonance (NMR)

CHEMICAL CHARACTERIZATION

Contact me at : jason.fauquet@umons.ac.be

Acknowledgement

We express our sincere gratitude to Prof. R. Wattiez from the Proteomics and Microbiology (ProMic) Unit for granting us access to the UMONS The Bioprofiling platform used for proteomic analysis was supported by the European Regional Development Fund and the Walloon Region, Belgium.

Poster created with canva.com